Revised: 29th Jul 2020 (3rd version)

- DIAGNOSTIC AGENT for Helicobacter Pylori Infection-

UBIT® Tablet 100 mg

< Urea (13C) Tablet >

Prescription-only Medication

Caution: Use only as directed by a physician.

Storage
Store below 30°C. Use immediately after opening SP package

Expiration date

Three years after the date of manufacture.

(The expiration date is indicated on the package.)

Singapore Registration no.	SIN13253P

COMPOSITION AND DESCRIPTION

1. Composition

Each tablet of UBIT Tablet 100mg contains 100 mg of urea (13C) and the following inactive ingredients: lactose hydrate, microcrystalline cellulose, corn starch, magnesium stearate, hypromellose, macrogol 6000, talc and titanium oxide.

2. Product description

UBIT Tablet 100mg is a white film-coated tablet.

Appearance	Diameter (mm)	Thickness (mm)	Weight (mg)	Code
(G) (C) (C) (C) (C) (C) (C) (C) (C) (C) (C	8.1	3.7	Approx.	OG73

INDICATIONS

In vivo diagnosis of Helicobacter pylori (H. pylori) infection

DOSAGE AND ADMINISTRATION

Usual adult dosage is one tablet [100 mg of urea (13C)] taken orally in a fasting condition.

<Recommended Procedure for 13C-Urea Breath Test>

- 1. Collect pre-dose breath sample.
- Instruct patient in a fasting condition to promptly (within 5 sec) swallow one tablet of UBIT Tablet 100mg with 100 mL of water without crushing or chewing the tablet.
- Instruct patient to remain in left lateral decubitus position for 5 min, and then to return to sitting position until collection of post-dose breath sample.
- Collect post-dose breath sample at 20 min after administration.
- Measure ¹³CO₂ level (ratio of ¹³CO₂ to ¹²CO₂) in each breath sample and calculate Δ¹³CO₂ (difference between pre-dose and post-dose ¹³CO₂ levels) for diagnosis of H. pylori infection.

PRECAUTIONS

1. Adverse Reactions

Adverse reactions were reported in 8 (0.7%) of a total of 1,144 patients who received UBIT Tablet 100mg or UBIT Granules 100mg (data available at the time of approval).

Adverse reactions were reported in 5 (0.14%) of 3,500 patients who received UBIT Granules 100mg in a post-marketing drug use results survey (data available at the time of completion of reasonination)

Signs and Symptoms/Incidence	<0.5% Note 1	Incidence unknown Note 2
Hypersensitivity Note 3		Rash, urticaria
Gastrointestinal	Abdominal distention, diarrhea, epigastric discomfort	Vomiting, nausea
Other	Increase in serum potassium level	

Note 1: Clinical-Trial Data

Note 2: Post-Authorization Data
The frequencies of reactions observed in the post-marketing period are considered unknown Note 3: If symptoms of hypersensitivity occur, the drug should be discontinued immediately

2. Use during Pregnancy, Delivery, or Lactation

The drug should be administered to pregnant, possibly pregnant, or lactating women only if the anticipated diagnostic benefit is thought to outweigh any potential risk. (The safety of this drug in pregnant and lactating women has not been established.)

3. Pediatric Use

The safety of this drug in low birth weight babies, newborns, suckling infants, infants, and children has not been established. (There is no clinical experience in low birth weight babies, newborns, suckling infants, and infants. There is little clinical experience in children.)

4. Interactions

Any treatment interfering with the H. pylori status or the urease activity can influence the urea breath test

5. Contraindications

Hypersensitivity to active substance or to any of the excipients

6. Overdose and Treatment

Due to the fact that 1 film-coated tablet contains only 100 mg of ¹³C-urea that reflects only a fraction of the natural daily food intake, an overdose is not to be expected.

<Pre><Precautions Concerning Use>

At time of ingestion, UBIT Tablet 100mg must be swallowed promptly (within 5 sec) without being crushed or chewed.

If the tablet disintegrates in the mouth, the diagnostic test result may be influenced by oral bacteria having urease activity, possibly causing a false-positive test result. When UBIT Tablet 100mg was placed in water, the tablet film coating began to flake in 5-8 sec.

<Pre><Precautions Concerning Diagnosis>

- (1) Diagnostic criteria: A Δ^{13} C value of 2.5% or higher at 20 min after administration of one tablet of **UBIT Tablet 100mg** [100 mg of urea (¹³C)] is considered a positive result for H. pylori infection.
- (2) Δ^{13} C values for breath samples in the 13 C-urea breath test should be determined by mass spectrometry or other appropriate methods with a comparable performance (e.g. infrared spectroscopy). It should be noted that a low CO₂ concentration (<1%) in the expired breath sample may cause poor reproducibility of the measured Δ^{13} C values¹⁾, possibly affecting the diagnostic test result at low Δ^{13} C values.



- (3) Practical considerations for the diagnosis of H. pylori infection: A falsenegative test result may be obtained if the ¹³C-urea breath test is performed during or immediately after discontinuation of therapy with medications that have bacteriostatic activity against H. pylori, such as proton pump inhibitors (omeprazole, lansoprazole, and sodium rabeprazole), antibiotics (amoxicillin hydrate, clarithromycin, and tetracycline), metronidazole, preparations, and ecabet sodium hydrate preparations having an inhibiting action on gastric urease activity. The urea breath tests should not be performed until 4 weeks without systemic antibacterial therapy and two (2) weeks after last dose of acid antisecretory agents. This is especially important after eradication therapy.
- (4) Diagnosis of H. pylori infection after eradication treatment (assessment of eradication therapy): Assessment of eradication therapy, when necessary, should be performed at least 4 weeks after completion of the therapy.
- (5) False-negative test results were obtained in approx. 2.3% (3 of 130 H. pyloriinfected subjects) in a Japanese Phase 3 study of UBIT Tablet 100mg. Results of other test methods should be used as a reference if patients show a negative test result despite clinical signs or symptoms suggestive of H. pylori infection.
- (6) Since this drug is a tablet, breath reaction could be delayed by various factors such as passage of the tablet through the stomach without being retained long enough for dissolution, possibly causing a false-negative test result2). Therefore, administration of UBIT Tablet 100mg to patients in whom a falsenegative test result is likely to occur, such as those who have undergone gastrectomy, should be performed carefully.
- (7) A false-positive test result may be obtained in patients with achlorhydria³⁾ or patients who are infected with gastric spiral microorganisms other than H. pylori that exhibit urease activity, such as Helicobacter heilmanii^{4,5)}.
- (8) No correlation has yet been established between the bacterial count of H. *pylori* and the 13 C-urea breath test result (the measured Δ^{13} C value).

PHARMACODYNAMICS/PHARMACOKINETICS **PHARMACOKINETICS**

ATC Code: VO4CX05

Serum concentration of exogenous urea (13C) following single oral administration of 100 mg of urea (13C) (one sachet of UBIT® Granules 100mg) to 5 healthy adult male volunteers in a fasting condition⁶⁾

H. pylori	T_{max}	C_{max}	t _{1/2}	AUC _{0-24hr}
antibody titer	(hr)	$(\mu g/mL)$	(hr)	$(\mu g \cdot hr/mL)$
Positive	0.6±0.2	2.3±0.5	7.7±1.7	21.3±5.4
Negative	1.1±1.1	3.0±0.8	5.7±1.6	24.1±13.4

Mean±S.D.

Data represent exogenous urea (13C) concentrations determined by subtracting the intrinsic urea (13C) concentration from the measured concentration.

Absorption

Almost all of the administered dose was absorbed in rats⁷⁾.

Plasma Protein Binding

6.6% in dogs8)

Major Metabolites and Major Metabolic Pathway

Urea (13 C) is not metabolized in the body^{6,7)}.

Excretion Route and Excretion Rate

Excretion route: Renal excretion⁶⁾

Excretion rate: When 14C-urea (50 mg as urea) dissolved in 20 mL water was administered once orally to 7 non-infected subjects after overnight fasting, the urinary excretion of ¹⁴C for 3 days post-dosing accounted for approx. 90% of the administered dose9).

CLINICAL STUDIES

In a Phase 3 study of UBIT Tablet 100mg, ¹³C-urea breath test was performed in 130 H. pylori-infected patients and 124 non-infected patients diagnosed by a combination of biopsy-based methods (FDA H. pylori Infection Assessment Criteria $1995)^{10)}$

Based on measurement of ¹³CO₂ levels in the expired breath by mass spectrometry employing a cut-off value of 2.5% at 20 min post-dosing, ¹³C-urea breath test using UBIT Tablet 100mg showed sensitivity of 97.7%, specificity of 98.4%, and accuracy of 98.0%.

Comparison of Δ^{13} C values measured by infrared spectroscopy with those measured by mass spectrometry at 2542 points indicated good correlation between the two measurement methods.

PHARMACOLOGY

H. pylori exhibits high urease activity.

The urea labeled with ¹³C as a tracer in UBIT Tablet 100mg is degraded to ¹³CO₂ and NH₃ by the urease activity of gastric *H. pylori* after ingestion, and the ¹³CO₂ that is formed is taken up by the blood and expired in the breath.

Based on the above principle, the 13 C-urea breath test is a method of diagnosing H. pylori infection by detecting the ¹³CO₂ in the expired breath and determining the Δ^{13} CO₂ after administration of urea (13 C).

PHYSICOCHEMISTRY

Nonproprietary name

Urea (13C) (JAN)

Chemical name

13C-urea

Structural formula:



Molecular formula

 $^{13}CH_4N_2O$

Molecular weight

61.05

Description

Urea (13C) occurs as colorless to white crystals or a crystalline powder and is odorless. It is highly soluble in water, soluble in ethanol (95), slightly soluble in acetonitrile, and practically insoluble in diethyl ether.

An aqueous solution (1 in 100) of Urea (13C) is neutral.

Melting point

133.5 - 135.5°C

PACKAGING

UBIT Tablet 100 mg: 1 x 10 tablets

REFERENCES

- 1) Kato M et al.: J. Jpn. So. Gastroenterol., 98 (7), 853, 2001.
- 2) Kato M et al.: Relationship between Gastric Mucosal Lesion and Helicobacter pylori - Latest Study Results -: Sentan-Igaku-Sha, Tokyo, 151-158, 1995.
- 3) Borriello SP et al.: J. Clin. Pathol., 38 (8), 946-953, 1985.
- 4) McNulty CAM. et al.: J. Clin. Pathol., 42 (6), 585-591, 1989.
- 5) Vaira D et al.: J. Clin. Pathol., 41 (7), 812-813, 1988.
- 6) Asaki S et al.: Jpn. J. Clin. Pharmacol. Ther., 27(3), 607-616, 1996.
- 7) Osumi T et al.: Jpn. Pharm. Ther., 24 (6), 1365-1374, 1996.
- 8) Watanabe J et al.: J. Pharm. Dyn., 7, 294-303, 1984.
- 9) Munster DJ et al.: Scand. J. Gastroenterol., 28 (8), 661-666, 1993.
- 10) FDA Center of Drug Evaluation and Research: Division of anti-infective drug products: Points to consider; Clinical development and labeling of anti-infective drug products, March 1995. Addendum (draft) Helicobacter pylori-associated peptic ulcer disease Indication # 25

Manufactured by:

Otsuka Pharmaceutical Co., Ltd.

Tokushima Factory

463-10, Kagasuno, Kawauchi-cho Tokushima-shi, Tokushima, Japan

Product Registrant:

Nagase Singapore (Pte) Ltd.

600 North Bridge Road, #11-01 Parkview Square,

Singapore 188778.

TEL: 65-6398-0088

FAX: 65-6398-0237